# Import and Export of Spectra Files

Vignette for the R package hyperSpec

Claudia Beleites <chemometrie@beleites.de> CENMAT and DI3, University of Trieste Spectroscopy · Imaging, IPHT Jena e.V.

March 3, 2015

# **Supported File Formats**

hyperSpec supports a number of file formats relevant for different types of spectroscopy. This is naturally only a subset of the file formats produced by different spectroscopic equipment. If you use hyperSpec with data formats not mentioned in this document, please send an email to Claudia Beleites <chemometrie@beleites.de>, so that this document can be updated. The information should include

- The type of spectroscopy
- Spectrometer model, manufacturer, and software
- The "native" file format (including a sample file)
- Description of relevant procedures to convert the file
- R code to import the data together with an example file that can actually be read by R.
- Documentation, particularly the description of the data format

If you need help finding out how to import your data, please search and eventually ask on Stack-exchange with tags [r] and [spectroscopy]. While I requiarly check these tags, consder dropping ma Claudia Beleites <chemometrie@beleites.de>an email in addition.

# Reproducing the Examples in this Vignette

The source code of this vignette including the spectra files are available as .zip file at hyperSpec's home page: http://hyperspec.r-forge.r-project.org/blob/fileio.zip
Note that some definitions are in file vignettes.defs.

# **Contents**

1. Introduction 2

2.	Creating a hyperSpec object with new 2.1. Creating a hyperSpec Object from a Data Matrix (Spectra Matrix)	3 3			
3. Reading Multiple files into one hyperSpec object					
4.	ASCII files 4.1. ASCII files with samples in columns 4.2. JCAMP-DX 4.3. Basic Atomic Spectra from NIST Tables 4.4. ASCII Export 4.5. ASCII Export 4.6. ASCII Export 4.7. ASCII Export 4.8. ASCII Export 4.9. ASCII Export 4.9. ASCII Export 4.0. ASCII Export 4.0. ASCII Export 4.0. ASCII Export	8			
5.	Binary file formats         5.1. Matlab Files       5.1.1. Matlab Export         5.1.2. Import of Matlab files written by Cytospec         5.2. ENVI Files       5.2.1. ENVI Export         5.3. spc Files       5.3. spc Files	8 9 9 10 10 10			
6.	Manufacturer-Specific Discussion of File Import  6.1. Manufacturer Specific Import Functions  6.2. Bruker FT-IR Imaging  6.3. Nicolet FT-IR Imaging  6.4. Varian/Agilent FT-IR Imaging  6.5. Kaiser Optical Systems Raman  6.5.1. Kaiser Optical Systems ASCII Files  6.5.2. Kaiser Optical Systems Raman Maps  6.6. Renishaw Raman  6.6.1. Renishaw ASCII data  6.7. Horiba / Jobin Yvon (e.g. LabRAM)  6.8. Witec	14 14 14 15 15 15 16 16 17			
7.	Writing your own Import Function 7.1. A new ASCII Import Function: scan.txt.PerkinElmer	18 19 20 22			
Α.	File Import Functions by Format	23			
В.	File Import Functions by Manufacturer	24			
C.	File Import Functions by Spectroscopy	25			

# 1. Introduction

This document describes how spectra can be imported into hyperSpec objects. Some possibilities to export hyperSpec objects as files are mentioned, too.

The most basic funtion to create *hyperSpec* objects is **new** ("hyperSpec") (section 2). It makes a *hyperSpec* object from data already in R's workspace. Thus, once the spectra are imported into R, conversion to *hyperSpec* objects is straightforward.

In addition, hyperSpec comes with predifined import functions for different data formats. This document divides the discussion into dealing with ASCII files (section 4, p. 5) and binary file formats (section 5, p. 8). If data export for the respective format is possible, it is discussed in the same sections. As sometimes the actual data written by the spectrometer software exhibits peculiarities, hyperSpec offers several specialized import functions. These are in general named after the data format followed by the manufacturer (e. g. read.ENVI.Nicolet).

Overview lists of the directly supported file formats are in the appendix: sorted by file format (appendix A, p. 23), manufacturer (appendix B, p. 24), and by spectroscopy (appendix C, p. 25).

# 2. Creating a hyperSpec object with new

To create a *hyperSpec* object from data in R's workspace, use:

```
> spc <- new ("hyperSpec", spc, wavelength, data, labels)
```

With the arguments:

spc the spectra matrix (may also be given as matrix inside column \$spc of data)

wavelength the wavelength axis vector

data the extra data (possibly already including the spectra matrix in column spc)

a list with the proper labels. Do not forget the wavelength axis label in \$.wavelength

and the spectral intensity axis label in \$spc.

Thus, once your data is in R's workspace, creating a hyperSpec object is easy. I suggest wrapping the code to import your data and the line joining it into a hyperSpec object by your own import function. You are more than welcome to contribute such import code to hyperSpec. Secion 7, (p. 18) discusses examples of custom import functions.

## 2.1. Creating a hyperSpec Object from a Data Matrix (Spectra Matrix)

As spectra matices are the internal format of hyperSpec, the consructor can directly be used:

```
> spc <- new ("hyperSpec", spc, wavelength, data, labels)
```

# 2.2. Creating a hyperSpec Object from a Data Cube (Spectra Array)

Roberto Moscetti asked how to convert a hyperspectral data cube into a hyperSpec object:

```
The problem is that I have a hypercube with the following dimensions: 67 \times 41 \times 256 \ y = 67 \ x = 41 \ wavelengths = 256
```

I do not know the way to import the hypercube.

Data cubes (i.e. 3-dimensional arrays of spectral data) result from spectal imaging measurements, where spectra are supplied for each pixel of an  $px.x \times px.y$  imaging area. They have 3 directions, usually x, y, and the spectral dimension.

The solution is to convert the array into a spectra matrix and have separate x and y coordinates.

Assume data is the data cube, and x, y and wl hold vectors with the proper x and y coordinates and the wavelengths:

```
> data <- array (1 : 24, 4 : 2)
> w1 <- c (550, 630)
> x <- c (1000, 1200, 1400)
> y <- c (1800, 1600, 1400, 1200)
> data
, , 1
     [,1] [,2] [,3]
[1,]
            6 10
[2.]
          7
              11
[3,]
       3
[4,]
           8
, , 2
     [,1] [,2] [,3]
[1,]
      14
           18
                22
[2,]
[3,]
      15
           19
                23
[4,]
      16
```

Such data can be converted into a *hyperSpec* object by:

```
> d <- dim (data)
> dim (data) <- c (d [1] * d [2], d [3])
> x <- rep (x, each = d [1])
> y <- rep (y, d [2])
> spectra <- new ("hyperSpec", spc = data,
+ data = data.frame (x, y), wavelength = wl)</pre>
```

If no proper coordinates (vectors x, y and w1) are available, they can be left out. In the case of x and y, map plotting will then be impossible, missing wavelengths will be replaced by column indices counting from 1 to d [3] automatically. Of course, such sequences (the row/column/pixel numbers) can be used instead of the original x and y as well:

```
> y <- seq_len (d [1])
> x <- seq_len (d [2])</pre>
```

Data cubes often come from spectral imaging systems that use an "image" coordinate system counting y from top to bottom. Note that this should accounted for in the decreasing order of the original y vector.

# 3. Reading Multiple files into one hyperSpec object

Many of the function described below will work on one file, even though derived functions such as read.spc.KaiserMap (see section 6.5.2, p. 15) may take care of measurements consisting of multiple files.

Usually, the most convenient way to import multiple files into one *hyperSpec* object is reading all files into a list of *hyperSpec* objects, and then collapseing this list into a single *hyperSpec* object:

```
> files <- Sys.glob ("spc.Kaisermap/*.spc")
> files <- files [seq (1, length (files), by = 2)] # import low wavenumber region only
> spc <- lapply (files, read.spc)
> length (spc)

[1] 54
> spc [[1]]
```

```
hyperSpec object
   1 spectra
   3 data columns
   1340 data points / spectrum
wavelength: x/"a. u." [numeric] 1 2 ... 1340
data: (1 rows x 3 columns)
   1. z: x/"a. u." [numeric] 1
   2. z.end: x/"a. u." [numeric] 1
   3. spc: Counts [matrix1340] 2782.7 2229.8 ... 932.02
> spc <- collapse (spc)
> spc
hyperSpec object
   54 spectra
   3 data columns
   1340 data points / spectrum
wavelength: x/"a. u." [numeric] 1 2 ... 1340
data: (54 rows x 3 columns)
   1. z: x/"a. u." [numeric] 1 1 ... 1 2. z.end: x/"a. u." [numeric] 1 1 ... 1
   3. spc: Counts [matrix1340] 2782.7 2678.6 ... 789.49
```

Note that in this particular case, the spectra are more efficiently read by read.spc.KaiserMap (see section 6.5.2, p. 15).

If you regularly import huge maps or images, writing a customized import function is highly encouraged. You may gain speed and memory by using the internal workhorse functions for the file import. In that case, please contact the package maintainer (Claudia Beleites <chemometrie@beleites.de>) for advise (contributions to hyperSpec are welcome and all authors are listed appropriately in the function help page's author section).

#### 4. ASCII files

Currently, hyperSpec provides two functions for general ASCII data import:

```
read.txt.long imports long format ASCII files, i.e. one intensity value per row
```

read.txt.wide imports wide format ASCII files, i.e. one spectrum per row

The import functions immediately return a hyperSpec object.

Internally, they use read.table, a very powerful ASCII import function. R supplies another ASCII import function, scan. scan imports numeric data matrices and is faster than read.table, but cannot import column names. If your data does not contain a header or it is not important and can safely be skipped, you may want to import your data using scan.

Note that R allows to use a variety of compressed file formats directly as ASCII files (for example, see section 6.6.1 on p. 16). Also, both read.txt.long and read.txt.wide accept connections instead of file names.

### 4.1. ASCII files with samples in columns

Richard Pena asked about importing another ASCII file type:

Triazine5\_31.txt file corresponds to X ray powder diffraction data (Bruker AXS). The native files data ".raw" are read with EVA software then they are converted into .uxd file with the File Exchange software (Bruker AXS). The .uxd file are opened with Excel

software and saved as .txt file, csv file (ChemoSpec) or xls.

The first and following columns corresponds to the angle diffraction and the intensity values of samples respectively.

This file thus differs from the ASCII formats discussed above in that the samples are actually in columns whereas *hyperSpec* expects them to be in rows. The header line gives the name of the sample. Import is straightforward, just the spectra matrix needs to be transposed to make a *hyperSpec* object:

```
> file <- read.table ("txt.t/Triazine 5_31.txt", header = TRUE, dec = ",", sep = "\t")
> triazine <- new ("hyperSpec", wavelength = file [,1], spc = t (file [, -1]),
                    data = data.frame (sample = colnames (file [, -1])),
                    labels = list (.wavelength = expression (2 * theta / degree),
                                   spc = "I / a.u."))
> triazine
hyperSpec object
   25 spectra
   2 data columns
   1759 data points / spectrum
wavelength: 2 * theta/degree [numeric] 5.0025 5.0173 ... 31.004
data: (25 rows x 2 columns)
   1. sample: [factor] DIV1208200 DIV1208300 ... VCA0106703
   2. spc: I / a.u. [matrix1759] 92 96 ... 163
> plot (triazine [1])
    0009
              10 12 14
                                   24 26 28 30
            8
                       16
                          18
                             20
                                22
         6
                          2θ/°
```

Witec also saves ASCII data with spectra in columns (Export  $\rightarrow$  Table), see 6.8.

## 4.2. JCAMP-DX

Limited import of JCAMP-DX files v. 4.24 [1] is available in function read.jdx. These files can contain multiple spectra, supported data formats are tabular XY..XY and X++(Y..Y).

```
> read.jdx ("jcamp-dx/shimadzu.jdx", encoding = "latin1", keys.hdr2data=TRUE)
```

```
hyperSpec object
    44 spectra
    11 data columns
    3401 data points / spectrum
wavelength: [numeric] 60.05 61.05 ... 548.3
data: (44 rows x 11 columns)
    1. spc: [matrix3401] 3294 NA ... 90 + NA
    2. file: [factor] jcamp-dx/shimadzu.jdx jcamp-dx/shimadzu.jdx ... jcamp-dx/shimadzu.jdx
    3. title: [numeric] 1 2 ... 44
    4. mw: [numeric] 261 249 ... 382
    5. molform: [factor] C11 H27 N O2 Si2 C9 H23 N O3 Si2 ... C25 H50 O2
    6. casregistryno: [factor] 72 - 18 - 4 56 - 45 - 1 ... 2442 - 49 - 1
    7. datatype: [factor] Mass Spectrum Mass Spectrum
    8. sampledescription: [factor] ...\11-12-12\pAS+OS Lauf1.qgd"\n20000 Da/s, ET 30ms, m/z 60-550\nPeak-Apex-Spektrum, H:
```

```
9. casname: [factor] L-Valin N,O-TMS2
                                               L-Serin 0,0'-TMS2 (X) ... Methyltetracosanoat (LRI C24)
   10. $retentionindex: [numeric] 884 925 ... 2400
   11. .format: [factor] (XY..XY) (XY..XY) ... (XY..XY)
> read.jdx ("jcamp-dx/virgilio.jdx")
hyperSpec object
   1 spectra
   1 data columns
   1216 data points / spectrum
wavelength: tilde(nu)/cm^-1 [numeric] 3030 3028 ... 600
data: (1 rows x 1 columns)
   1. spc: A [matrix1216] -3.1244e-05 0.0000e+00 ... 0
The last file has a slight inconsistenty between its meta data and spectroscopic data, causing a
message. However, the difference is minute compared to the intensities. If this is known in advance,
an appropriate tolerance can be chosen:
> read.jdx ("jcamp-dx/virgilio.jdx", ytol = 1e-9)
hyperSpec object
   1 spectra
   1 data columns
   1216 data points / spectrum
wavelength: tilde(nu)/cm^-1 [numeric] 3030 3028 ... 600
data: (1 rows x 1 columns)
   1. spc: A [matrix1216] -3.1244e-05 0.0000e+00 ... 0
To test the import, read a number of test files:
jcamp-dx/AMA1.DX
jcamp-dx/AMA2.DX
jcamp-dx/AMA3.DX
jcamp-dx/br_154_1.DX
jcamp-dx/BRUKAFFN.DX
jcamp-dx/BRUKPAC.DX
jcamp-dx/ISAS_MS1.DX
jcamp-dx/LABCALC.DX
jcamp-dx/PE1800.DX
jcamp-dx/testjose.dx
jcamp-dx/PE-IR/br_1.DX
jcamp-dx/PE-IR/br_2.DX
jcamp-dx/PE-IR/br_3.DX
jcamp-dx/PE-IR/br_4.DX
jcamp-dx/PE-IR/br_5.DX
jcamp-dx/PE-IR/fort_1.DX
jcamp-dx/PE-IR/fort_2.DX
jcamp-dx/PE-IR/fort_3.DX
jcamp-dx/PE-IR/fort_4.DX
jcamp-dx/PE-IR/fort_5.DX
jcamp-dx/PE-IR/lp_1.DX
jcamp-dx/PE-IR/1p_2.DX
jcamp-dx/PE-IR/lp_3.DX
jcamp-dx/PE-IR/1p_4.DX
jcamp-dx/PE-IR/lp_5.DX
jcamp-dx/GMD_20111121_MDN35_ALK_JCAMP.txt
```

#### Note

read.jdx.Shimadzu is deprecated.

## 4.3. Basic Atomic Spectra from NIST Tables

The NIST (National Institute of Standards and Technology) has published a data base of basic atomic emission spectra (see http://physics.nist.gov/PhysRefData/Handbook/periodictable.htm)[?] with emission lines tabulated in ASCII (HTML) files.

Here's an example how to extract the data of the Hg strong lines file:

```
> file <- readLines("NIST/mercurytable2.htm")</pre>
> #file <- readLines("http://physics.nist.gov/PhysRefData/Handbook/Tables/mercurytable2.htm")
> file <- file [- (1 : grep ("Intensity.*Wavelength", file) - 1)]
> file <- file [1 : (grep ("</pre>", file) [1] - 1)]
> file <- gsub ("<[^>]*>", "", file)
> file <- file [! grepl ("^[[:space:]]+$", file)]</pre>
> colnames <- file [1]
> colnames <- gsub ("[[:space:]][[:space:]]+", "\t", file [1])
> colnames <- strsplit (colnames, "\t")[[1]]</pre>
> if (colnames != c ("Intensity", "Wavelength (Å)", "Spectrum", "Ref. "))
    stop ("file format changed!")
> tablestart <- grep ("^[[:blank:]]*[[:alpha:]]+$", file) + 1</pre>
            <- c (tablestart [-1] - 2, length (file))
> tables <- list ()
> for (t in seq_along (tablestart)){
    tmp <- file [tablestart [t] : tableend [t]]</pre>
    tables [[t]] <- read.fwf (textConnection (tmp), c (5, 8, 12, 15, 9))
    \verb|colnames| (tables [[t]]) <- c("Intensity", "persistent", "Wavelength", "Spectrum", "Ref. ")| \\
    tables [[t]]$type <- gsub ("[[:space:]]", "", file [tablestart [t] - 1])</pre>
+ }
> tables <- do.call (rbind, tables)
> levels (tables$Spectrum) <- gsub (" ", "", levels (tables$Spectrum))
> Hg.AES <- list ()
> for (s in levels (tables$Spectrum))
    Hg.AES [[s]] <- new ("hyperSpec", wavelength = tables$Wavelength [tables$Spectrum == s],</pre>
                          spc = tables$Intensity [tables$Spectrum == s],
                          data = data.frame (Spectrum = s),
                          label = list (.wavelength = expression (lambda / ring (A)),
                                         spc = "I"))
> plot (collapse (Hg.AES), lines.args = list (type = "h"), col = 1 : 2)
```

# 4.4. ASCII Export

ASCII export can be done in wide and long format using write.txt.long and write.txt.wide. If you need a specific header or footer, use R's functions for writing files: write.table, write, cat and so on offer fine-grained control of writing ASCII files.

# 5. Binary file formats

#### 5.1. Matlab Files

Matlab files can be read and written using the package R.matlab[2], which is available at CRAN and can be installed by install.packages ("R.matlab").

```
spc.mat <- readMat ("spectra.mat")</pre>
```

If the .mat file was saved with compression, the additional package *Rcompression* is needed. It can be installed from omegahat:

```
install.packages("Rcompression", repos = "http://www.omegahat.org/R")
```

See the documentation of R.matlab for more details and possibly needed further packages.

readMat imports the .mat file's contents as a list. The variables in the .mat file are properly named elements of the list. The *hyperSpec* object can be created using new, see 2 (p. 3).

Again, you probably want to wrap the import of your matlab files into a function.

## 5.1.1. Matlab Export

R.matlab's function writeMat can be used to write R objects into .mat files. To save an hyperSpec object x for use in Matlab, you most likely want to save:

- the wavelength axis as obtained by wl (x),
- the spectra matrix as obtained by x [[]], and
- possibly also the extra data as obtained by x\$..
- as well as the axis labels labels (x).
- Alternatively, x\$. yields the extra data together with the spectra matrix.

However, it may be convenient to transform the saved data according to how it is needed in Matlab. The functions as.long.df and as.wide.df may prove useful for reshaping the data.

#### 5.1.2. Import of Matlab files written by Cytospec

A custom import function for .mat files written by Cytospec is available.

Note that Cytospec files can contain multiple versions of the data, the so-called blocks. The block to be read can be specified with the block argument. TRUE will read all blocks into a list:

```
> read.cytomat ("mat.cytospec/cytospec.mat", blocks = TRUE)
```

```
hyperSpec object
   55 spectra
   5 data columns
   981 data points / spectrum
wavelength: [numeric] 499.12 501.77 ... 3100
data: (55 rows x 5 columns)
   1. x: [integer] 4 5 ... 7
   2. y: [integer] 1 1 ... 11
   3. file: [factor] mat.cytospec/cytospec.mat mat.cytospec/cytospec.mat ... mat.cytospec/cytospec.mat
   4. block: [integer] 1 1 ... 1
   5. spc: [matrix981] 2112.9 2114.3 ... 2323.3
[[2]]
hyperSpec object
   55 spectra
   5 data columns
   981 data points / spectrum
wavelength: [numeric] 499.12 501.77 ... 3100
data: (55 rows x 5 columns)
   1. x: [integer] 4 5 ... 7
   2. y: [integer] 1 1 ... 11
   \textbf{3. file:} \quad \texttt{[factor]} \quad \texttt{mat.cytospec/cytospec.mat} \quad \texttt{mat.cytospec/cytospec.mat} \quad \dots \quad \texttt{mat.cytospec/cytospec.mat} \\
   4. block: [integer] 2 2 ... 2
   5. spc: [matrix981] 58.472 59.024 ... 262.5
```

#### 5.2. ENVI Files

ENVI files are binary data accompanied by an ASCII header file. *hyperSpec*'s function read.ENVI can be used to import them. Usually, the header file name is the same as the binary data file name with the suffix replaced by .hdr. Otherwise, the header file name can be given via parameter *headerfile*.

As we experienced missing header files (Bruker's Opus software frequently produced header files without any content), the data that would usually be read from the header file can also be handed to read. ENVI as a list in parameter *header*. Arguments given in *header* replace corresponding entries of the header file. The help page gives details on what elements the list should contain, see also the discussion of ENVI files written by Bruker's OPUS software (section 6.2, p. 14).

Here is how to use read. ENVI:

```
> spc <- read.ENVI ("ENVI/example2.img")
> spc

hyperSpec object
    0 spectra
    3 data columns
    1738 data points / spectrum
wavelength: [numeric] 649.90 651.83 ... 3999.7
data: (0 rows x 3 columns)
    1. x: [integer]
    2. y: [integer]
    3. spc: [matrix1738]
```

Please see also the manufacturer specific notes in section 6.1, p. 14.

# 5.2.1. ENVI Export

Use package caTools or rgdal with GDAL for writing ENVI files.

#### 5.3. spc Files

Thermo Galactic's .spc file format[3] can be imported by read.spc.

A variety of sub-formats exists. *hyperSpec*'s import function read.spc does *not* support the old file format that was used before 1996. In addition, no test data with *w planes* was available — thus the import of such files could not be tested. If you come across such files, please contact the package maintainer (Claudia Beleites <chemometrie@beleites.de>).

Here are some tests using Thermo Galactic's example files:

```
> ## old format files stop with an error:
> old <- paste ("spc", c ('CONTOUR.SPC', 'DEMO 3D.SPC', 'LC DIODE ARRAY.SPC'), sep = "/")
> for (f in old)
+ try (read.spc (f))
> ## all other files should be good for import
> other <- setdiff (Sys.glob ("spc/*.[sS][pP][cC]"), old)
> for (f in other){
+ spc <- read.spc (f)
+
+ if (is (spc, "hyperSpec"))
+ cat (f, ": ", nrow (spc), " spectrum(a), ", nwl (spc), " data pts / spc.\n", sep = "")
+ else</pre>
```

```
cat (f, ": list of ", length (spc), " spectra, ",
            paste (range (sapply (spc, nwl)), collapse = " - "),
            " data pts / spc\n", sep = "")
+ }
spc/BARBITUATES.SPC: list of 286 spectra, 4 - 101 data pts / spc
spc/barbsvd.spc: list of 286 spectra, 4 - 101 data pts / spc
spc/BENZENE.SPC: 1 spectrum(a), 1842 data pts / spc.
spc/DRUG SAMPLE_PEAKS.SPC: list of 6 spectra, 80 - 253 data pts / spc
spc/DRUG SAMPLE.SPC: list of 400 spectra, 2 - 254 data pts / spc
spc/FID.SPC: 1 spectrum(a), 8192 data pts / spc.
spc/HCL.SPC: 1 spectrum(a), 8361 data pts / spc.
spc/HOLMIUM.SPC: 1 spectrum(a), 901 data pts / spc.
spc/IG_BKGND.SPC: 1 spectrum(a), 4096 data pts / spc.
spc/IG_MULTI.SPC: 10 spectrum(a), 4096 data pts / spc.
spc/IG_SAMP.SPC: 1 spectrum(a), 4645 data pts / spc.
spc/KKSAM.SPC: 1 spectrum(a), 751 data pts / spc.
spc/POLYR.SPC: 1 spectrum(a), 1844 data pts / spc.
spc/POLYS.SPC: 1 spectrum(a), 1844 data pts / spc.
spc/SINGLE POLYMER FILM.SPC: 1 spectrum(a), 1844 data pts / spc.
spc/SPECTRUM WITH BAD BASELINE.SPC: 1 spectrum(a), 1400 data pts / spc.
spc/TOLUENE.SPC: 1 spectrum(a), 801 data pts / spc.
spc/TriVista-linear.spc: 1 spectrum(a), 1149 data pts / spc.
spc/TriVista-normal.spc: 1 spectrum(a), 1024 data pts / spc.
spc/TUMIX.SPC: 1 spectrum(a), 1775 data pts / spc.
spc/TWO POLYMER FILMS.SPC: 1 spectrum(a), 1844 data pts / spc.
spc/Witec-timeseries.spc: 25 spectrum(a), 1024 data pts / spc.
spc/XYTRACE.SPC: 1 spectrum(a), 3469 data pts / spc.
```

The header and subheader blocks of spc files store additional information of pre-defined types (see the file format specification[3]). Further information can be stored in the so-called log block at the end of the file, and should be in a key-value format (although even the official example files do not always). This information is often useful (Kaiser's Hologram software e.g. stores the stage position in the log block).

read.spc has four arguments that allow fine-grained control of storing such information in the hyperSpec object:

keys.hdr2data parameters from the spc file and subfile headers that should become extra data columns

keys.log2data parameters from the spc file log block that should become extra data columns

keys.\*2log parameters are deprecated because the logbook itself is depecated

The value of these arguments can either be logical (amounting to either use all or none of the information in the file) or a character vector giving the names of the parameters that should be used. Note that the header file field names are always lowercase.

Here's how to find out what extra information could be read from the header and log:

> read.spc ("spc.Kaisermap/ebroAVII.spc", keys.hdr2data = TRUE)

```
hyperSpec object

1 spectra

33 data columns

1340 data points / spectrum

wavelength: x/"a. u." [numeric] 1 2 ... 1340

data: (1 rows x 33 columns)

1. z: x/"a. u." [numeric] 1

2. z.end: x/"a. u." [numeric] 1

3. ftflgs: [logical] FALSE

4. fexper: [factor] General

5. fexp: [integer] -128
```

```
6. fnpts: [integer] 1340
   7. ffirst: [numeric] 1
   8. flast: [numeric] 1340
   9. fnsub: [integer] 1
   10. fxtype: [character] x/"a. u."
   11. fytype: [character] Counts
   12. fztype: [character] x/"a. u."
13. fpost: [integer] 0
   14. fdate: [POSIXct, POSIXt] 1253590860
   15. fres: [character]
   16. fsource: [character]
   17. fspare: [numeric] 0
   18. fcmnt: [character]
   19. fcatxt: [character]
20. flogoff: [integer] 5904
   21. fmods: [integer] 0
   22. fprocs: [integer] 0
23. flevel: [integer] 0
   24. fsampin: [integer] 0
   25. ffactor: [numeric] 0
26. fmethod: [character]
   27. fzinc: [numeric] 0
   28. fwplanes: [integer] 0
   29. fwinc: [numeric] 0
   30. fwtype: [character] x/"a. u."
   31. .last.read: [numeric] 512 32. subfiledir: [numeric] 0
   33. spc: Counts [matrix1340] 2782.7 2229.8 ... 932.02
> read.spc ("spc.Kaisermap/ebroAVII.spc", keys.log2data = TRUE)
hyperSpec object
   1 spectra
   54 data columns
   1340 data points / spectrum
wavelength: x/"a. u." [numeric] 1 2 ... 1340
data: (1 rows x 54 columns)
   1. z: x/"a. u." [character] 1
   2. z.end: x/"a. u." [character] 1
3. Grams_File_Name: [character] d:\beleites\ebro\Map 20090921 180944\ebroAVII.spc
   4. HoloGRAMS_File_Name: [character] Unknown
   5. Acquisition_Date_Time: [character] 22.09.2009 03:41:30
   6. Lambda: [character] Low
   7. Accuracy_Mode: [character] High Speed
   8. Dark_subtracted: [character] Yes
9. Dark_File_Name: [character] ebroAVGC.drk
   10. Auto_New_Dark_Curve: [character] No
   11. Background_subtracted: [character] No
12. Background_File_Name: [character] <None>
   13. Intensity_Corrected: [character] Yes
   14. Intensity_Calibration_Available: [character] Yes
   15. Intensity_Correction_File: [character] c:\hologram\calibration\intensity\20090609aa.icl
   16. Intensity_Correction_Threshold: [character] 0,00%
   17. Intensity_Source_Correction: [character] No
   18. Intensity_Source_Correction_File: [character] <None>
   19. Comment: [character] <None>
   20. Cosmic_Ray_Filtering: [character] Yes
   21. Total_Cosmic_Count: [character] 38
   22. Exposure_Length: [character] 10000
   23. Accumulations: [character] 1
   24. Accumulation_Method: [character] Averaged
   25. Calibration_File: [character] C:\HoloGRAM\calibration\Wavelength\20090716ab.wcl
   26. Comment.1: [character]
   27. Temperature_Status: [character] At temperature
   28. Temperature: [character] -85,50
   29. HoloGRAMS_File_Version: [character] 4.1
```

```
31. Operator: [character] unknown
   32. Stage_X_Position: [character] 360
   33. Stage_Y_Position: [character] 1100 34. Stage_Z_Position: [character] -16
   35. AutoFocusUsed: [character] Nein
   36. WLInterval: [character] 0,12
37. CalInterval: [character] 0,90
   38. FFTFillFactor: [character] None
   39. FFTApT: [character] None
   40. SamplingMethod: [character] Linear
   41. Has_MultiPlex_Laser: [character] No
   42. External_Trigger: [character] No 43. Laser_Wavelength: [character] 785,21
   44. Default_Laser_Wavelength: [character] 785,21
   45. Laser_Tracking: [character] False
   46. Laser_Block_Active: [character] No
   47. Pixel_Fill_minimum: [character] 26
   48. Pixel_Fill_maximum: [character] 16762
   49. Binning_Start: [character] 24
   50. Binning_End: [character] 41
   51. NumPoints: [character] 1340
   52. First: [character] 1,000
   53. last: [character] 1340,000
   54. spc: Counts [matrix1340] 2782.7 2229.8 ... 932.02
.spc files may contain multiple spectra that do not share a common wavelength axis. In this case,
read.spc returns a list of hyperSpec objects with one spectrum each. collapse may be used to
combine this list into one hyperSpec object:
> barbiturates <- read.spc ("spc/BARBITUATES.SPC")</pre>
> save (barbiturates, file = "barbiturates.rda")
> class (barbiturates)
[1] "list"
> length (barbiturates)
[1] 286
> barbiturates <- do.call (collapse, barbiturates)
> barbiturates <- orderwl (barbiturates)
> barbiturates
hyperSpec object
   286 spectra
   3 data columns
   375 data points / spectrum
wavelength: frac(m, z)/frac(u, e) [numeric] 25.95 26.05 ... 244.05
data: (286 rows x 3 columns)
   1. z: t/min [numeric] 4.0272 4.0341 ... 5.9978
   2. z.end: t/min [numeric] 4.0272 4.0341 ... 5.9978
   3. spc: I/"a. u." [matrix375] NA NA ... NA + NA
> barbiturates [[1:10, , 25 ~ 30]]
      25.95 26.05 26.15 26.95 27.05 27.15 28.05 28.15 29.05 29.15 29.95
                     NA
                                 562
                                             NA 11511 6146
 [1,]
         NA
               NA
                           NA
                                       NA
 [2,]
               NA
                     NA
                           NA
                                 NA
                                       618 10151 NA 5040
                                                                       NA
         NΑ
                                                                 NΑ
 [3,]
         NA
               NA
                     NΑ
                           NΑ
                                 638
                                       NA
                                              NA 10722 5253
                                                                       NΑ
 [4,]
         NA
               NA
                     NA
                           NA
                                  NA
                                        NA 10548
                                                   NA
                                                        5865
                                                                 NA
                                                                       NA
 [5,]
         NΑ
               NA
                     NA
                           NΑ
                                  NA
                                        NA
                                             NA 10519
                                                        4664
                                                                 NΑ
                                                                       NΔ
 [6,]
               NA
                     NA
                           796
                                        NA 10519
                                                  NA 5110
         NA
                                  NA
 [7,]
         NA
               NA
                     NA
                           NA
                                  NA
                                        NA 10096
                                                    NA
                                                        4769
                                                                 NA
                                                                     907
 [8,]
         NA
               NA
                     NΑ
                           NΑ
                                  NA
                                       NA NA 10929
                                                       5400
                                                                 NA
                                                                      NΑ
 [9,]
         NA
               NA
                     NA
                           NA
                                 NA
                                       NA 10235
                                                 NA 4930
                                              NA 10663 4690
[10,]
                     NA
                           NA
                                 NA
                                       NA
                                                                     799
         NA
               NA
                                                                 NΑ
```

30. File\_Type: [character] Hol

**Deriving manufacturer specific import filters.** Please note that future changes inside the read.spc function are likely to occur. However, if you just post-process the *hyperSpec* object returned by read.spc, you should be fine.

## 6. Manufacturer-Specific Discussion of File Import

#### 6.1. Manufacturer Specific Import Functions

Many spectrometer manufacturers provide a function to export their spectra into ASCII files. The functions discussed above are written in a very general way, and are highly customizable. I recommend wrapping these calls with the appropriate settings for your spectra format in an import function. Please consider contributing such import filters to *hyperSpec*: send me the documented code (for details see the box at the beginning of this document). If you are able to import data of any format not mentioned in this document (even without the need of new converters), please let me know (details again in the box at the beginning of this document).

# 6.2. Bruker FT-IR Imaging

We use read.ENVI to import IR-Images collected with a Bruker Hyperion spectrometer with OPUS software. As mentioned above, the header files are frequently empty. We found the necessary information to be:

```
> header <- list (samples = 64 * no.images.in.row,
+ lines = 64 * no.images.in.column,
+ bands = no.data.points.per.spectrum,
+ 'data type' = 4,
+ interleave = "bip")</pre>
```

No spatial information is given in the ENVI header (if written). The lateral coordinates can be setup by specifying origin and pixel size for x and y directions. For details please see the help page.

The proprietary file format of the Opus software is not yet supported.

#### 6.3. Nicolet FT-IR Imaging

Also Nicolet saves imaging data in ENVI files. These files use some non-standard keywords in the header file that should allow to reconstruct the lateral coordinates as well as the wavelength axes and units for wavelength and intensity axis. *hyperSpec* has a specialized function read.ENVI.Nicolet that uses these header entries.

It seems that the position of the first spectrum is recorded in µm, while the pixel size is in mm. Thus a flag *nicolet.correction* is provided that divides the pixel size by 1000. Alternatively, the correct offset and pixel size values may be given as function arguments.

```
> spc <- read.ENVI.Nicolet ("ENVI/example2.img", nicolet.correction = TRUE)
> spc ## dummy sample with all intensities zero
hyperSpec object
    0 spectra
    3 data columns
    1738 data points / spectrum
wavelength: [numeric] 649.90 651.83 ... 3999.7
data: (0 rows x 3 columns)
    1. x: [numeric]
    2. y: [numeric]
    3. spc: [matrix1738]
```

# 6.4. Varian/Agilent FT-IR Imaging

Agilent (Varian) uses a variant of ENVI (with binary header). A specialized form of read.ENVI will be coming soon.

## 6.5. Kaiser Optical Systems Raman

Spectra obtained using Kaiser's Hologram software can be saved either in their own .hol format and imported into Matlab (from where the data may be written to a .mat file readable by R.matlab's readMat. Hologram can also write ASCII files and .spc files. We found working with .spc files the best option.

The spectra are usually interpolated by Hologram to an evenly spaced wavelength (or  $\Delta \tilde{\nu}$ ) axis unless the spectra are saved in a by-pixel manner. In this case, the full spectra consist of two files with consecutive file names: one for the low and one for the high wavenumber region. See the example for .spc import.

#### 6.5.1. Kaiser Optical Systems ASCII Files

The ASCII files are long format that can be imported by read.txt.long (see section 4, p. 5).

We experienced two different problems with these files:

- 1. If the instrument computer's locale is set so that also the decimal separator is a comma, commas are used both as decimal and as column separator.
- 2. Values with a decimal fraction of 0 are written with decimal separator but no further digits (e.g. 2,). This may be a problem for certain conversion functions (read.table works fine, though).

Thus care must be taken:

```
> ## 1. import as character
> tmp <- scan ("txt.Kaiser/test-lo-4.txt", what = rep ("character",4), sep = ",")
> tmp <- matrix (tmp, nrow = 4)
> ## 2. concatenate every two columns by a dot
> wl <- apply (tmp [1:2, ], 2, paste, collapse = '.')
> spc <- apply (tmp [3:4, ], 2, paste, collapse = '.')
> ## 3. convert to numeric and create hyperSpec object
> spc <- new ("hyperSpec", spc = as.numeric (spc), wavelength = as.numeric (wl))</pre>
```

#### 6.5.2. Kaiser Optical Systems Raman Maps

hyperSpec provides the function read.spc.KaiserMap to easily import spatial collections of .spc files written by Kaiser's Hologram software. The filenames of all .spc files to be read into one hyperSpec object can be provided either as a character vector or as a wildcard expression (e.g. "path/to/files/\*.spc").

The data for the following example was saved with wavelength axis being camera pixels rather than Raman shift. Thus two files for each spectrum were saved by Hologram. Thus, a file name pattern is difficult to give and a vector of file names is used instead:

```
> files <- Sys.glob ("spc.Kaisermap/*.spc")
> spc.low <- read.spc.KaiserMap (files [seq (1, length (files), by = 2)])
> spc.high <- read.spc.KaiserMap (files [seq (2, length (files), by = 2)])
> wl (spc.high) <- wl (spc.high) + 1340
> spc
```

```
hyperSpec object
  1 spectra
  1 data columns
  2110 data points / spectrum
wavelength: [numeric] 121.5 122.4 ... 2019.6
data: (1 rows x 1 columns)
  1. spc: [matrix2110] 1202.51 770.35 ... 141.01
```

#### 6.6. Renishaw Raman

Renishaw's Wire software comes with an file format converter. This program can produce a long ASCII format, .spc, or .jdx files.

We experienced that the conversion to .spc is *not* fully reliable: maps were saved as depth profile, loosing all spatial information. In addition, an evenly spaced wavelength axis was produced, although this was de-selected in the converter. We therefore recommend using the ASCII format. Otherwise the import using read.spc worked.

#### 6.6.1. Renishaw ASCII data

An optimized import function for the ASCII files is available: scan.txt.Renishaw. The file may be compressed via gzip, bzip2, xz or lzma. zip compressed files are read via scan.zip.Renishaw. The ASCII files can easily become very large, particularly with linefocus- or streamline imaging. scan.txt.Renishaw provides two mechanisms to avoid running out of memory during data import. The file may be imported in chunks of a given number of lines (see the last example). scan.txt.Renishaw can calculate the correct number of wavelengths (i.e. data points per spectrum) if the system command wc is available on your computer.

In addition, the processing of the long ASCII format into the spectra matrix is done by reshaping the vector of intensities into a matrix. This process does not allow any missing values in the data. Therefore it is not possible to import multi-spectra files with individually "zapped" spectra using scan.txt.Renishaw.

The second argument to scan.txt.Renishaw decides what type of experiment is imported. Supported types are:

```
"xyspc"
                  maps, images, multiple spectra with x and y coordinates (default)
"spc"
                  single spectrum
"depth", "zspc" depth series
"ts"
                  time series
Instead of a file name, scan.txt.Renishaw accepts also a connection.
> paracetamol <- scan.txt.Renishaw ("txt.Renishaw/paracetamol.txt", "spc")</pre>
> paracetamol
hyperSpec object
   1 spectra
   1 data columns
   4064 data points / spectrum
wavelength: Delta * tilde(nu)/cm^-1 [numeric] 96.787 98.143 ... 3200.1
data: (1 rows x 1 columns)
   1. spc: I / a.u. [matrix4064] 2056.5 2224.8 ... 299.23
> save (paracetamol, file = "paracetamol.rda")
> scan.txt.Renishaw ("txt.Renishaw/laser.txt.gz", "ts")
```

```
hyperSpec object
   84 spectra
   2 data columns
   140 data points / spectrum
wavelength: Delta * tilde(nu)/cm^-1 [numeric] -199.08 -196.90 ... 99.934
data: (84 rows x 2 columns)
   1. t: t / s [numeric] 0 2 ... 5722
   2. spc: I / a.u. [matrix140] 29.801 32.093 ... 81.3
Very large files can be read in chunks to save memory:
> scan.txt.Renishaw ("txt.Renishaw/chondro.txt", nlines = 1e5, nspc = 875)
hyperSpec object
  875 spectra
   3 data columns
   1272 data points / spectrum
wavelength: Delta * tilde(nu)/cm^-1 [numeric] 601.62 602.66 ... 1802.2
data: (875 rows x 3 columns)
   1. y: y/(mu * m) [numeric] -4.77 -4.77 ... 19.23
   2. x: x/(mu * m) [numeric] -11.55 -10.55 ... 22.45
   3. spc: I / a.u. [matrix1272] 501.72 518.53 ... 151.92 + NA
R accepts a variety of compressed file formats for ASCII files:
> scan.txt.Renishaw ("txt.Renishaw/chondro.gz")
> scan.txt.Renishaw ("txt.Renishaw/chondro.xz")
> scan.txt.Renishaw ("txt.Renishaw/chondro.lzma")
> scan.txt.Renishaw ("txt.Renishaw/chondro.gz")
> scan.txt.Renishaw ("txt.Renishaw/chondro.bz2")
> scan.zip.Renishaw ("txt.Renishaw/chondro.zip")
```

# 6.7. Horiba / Jobin Yvon (e.g. LabRAM)

Horiba's Labspec software (e.g. LabRAM spectrometers) saves spectra in a wide ASCII format which is read by read.txt.Horiba, e.g.:

Note that Labspec .txt files can contains lots of spectra with zero intensity: Labspec saves a complete rectangular grid even if only part of a map was measured. These spectra are by removed (remove.zerospc = TRUE).

For convenience, functions to further wrappers to import maps (read.txt.Horiba.xy) and time series (read.txt.Horiba.t) are provided.

#### 6.8. Witec

> read.spc ("spc/Witec-timeseries.spc")

The Witec project software supports exporting spectra as Thermo Galactic .spc files.

```
hyperSpec object
25 spectra
3 data columns
1024 data points / spectrum
wavelength: Delta * tilde(nu)/cm^-1 [numeric] 579.89 582.97 ... 3153.3
data: (25 rows x 3 columns)
1. z: T/(degree * C) [numeric] 0 0 ... 0
2. z.end: T/(degree * C) [numeric] 5.7509e-42 5.7509e-42 ... 5.7509e-42
3. spc: I/"a. u." [matrix1024] 1004 1008 ... 1006
```

.spc is in general the recommended format for hyperSpec import, but this export option is not supported for imaging data.

Imaging data can be exported as ASCII X and Y files (Save ASCII X and Save ASCII Y). These can be read by scan.dat.Witec:

```
> scan.dat.Witec (filex = "txt.Witec/WitecExample-x.dat", points.per.line = 10, lines.per.image = 10)
hyperSpec object
  100 spectra
  4 data columns
  1024 data points / spectrum
wavelength: [numeric] 106.39 110.76 ... 3377.2
data: (100 rows x 4 columns)
  1. spc: [matrix1024] 1494 1486 ... 952
  2. x: [integer] 1 2 ... 10
  3. y: [integer] -1 -1 ... -10
  4. filename: filename [character] txt.Witec/WitecExample-y.dat txt.Witec/WitecExample-y.dat ... txt.Witec/WitecExample-y.dat
```

Note that the Y data files also contain a wavelength information, but (at least Witec Project 2.10) this information is always wavelength in nm, not Raman shift in wavenumbers: this is provided by the X data file only.

Another option is Witec's txt table ASCII export (Export → Table), which produces ASCII files with each row corresponding to one wavelength. Such files can be read with scan.txt.Witec:

```
hyperSpec object
1 spectra
2 data columns
51200 data points / spectrum
wavelength: [numeric] 529.14 603.00 ... 1702
data: (1 rows x 2 columns)
1. spc: [matrix51200] 529.14 603.00 ... 1702
2. filename: filename [character] txt.Witec/Witec.txt
```

> scan.dat.Witec ("txt.Witec/Witec.txt")

scan.txt.Witec assumes 1024 wavelengths, but other values may be given as parameter nwl. NULL indicates that the number of wavelengths should be determined automatically.

#### 7. Writing your own Import Function

This section gives examples how to write import functions. The first example implements an import filter for an ASCII file format basically from scratch. The second example shows how to implement more details for an already existing import filter.

#### 7.1. A new ASCII Import Function: scan.txt.PerkinElmer

The raw spectra of the flu data set (see also the respective vignette) are in PerkinElmer's ASCII file format, one spectrum per file.

We need a function that automatically reads all files specified by a pattern, such as \*.txt. In order to gain speed, the spectra matrix should be preallocated after the first file is read.

A short examination of the files (flu\*.txt in directory txt.PerkinElmer) reveals that the actual spectrum starts at line 55, after a line containing #DATA. For now, no other information of the files is to be extracted. It is thus easier to skip the first 54 lines than searching for the line after #DATA.

A fully featured import function should support:

- Reading multiple files by giving a pattern
- hand further arguments to scan. This comes handy in case the function is used later to import other data types.
- Also skipping 54 lines would be a weird default, so we rather require it to be given explicitly.
- The same applies for the axis labels: they should default to reasonable settings for fluorescence spectra, but it should be possible to change them if needed.
- The usual log entry arguments should be supplied.
- A sanity check should be implemented: stop with an error if a file does not have the same wavelength axis as the others.
- Finally, if no file can be found, an empty *hyperSpec* object is a reasonable result: There is no need to stop with an error, but it is polite to issue an additional warning.

```
_ scan.txt.PerkinElmer.R
scan.txt.PerkinElmer <- function (files = "*.txt", ..., label = list ()) {</pre>
  ## set some defaults
 long <- list (files = files, ..., label = label)</pre>
  label <- modifyList (list (.wavelength = expression (lambda / nm),
                               spc = expression (I[f1] / "a.u.")),
                         label)
  ## find the files
  files <- Sys.glob (files)
  if (length (files) == 0){
    warning ("No files found.")
    return (new ("hyperSpec"))
  ## read the first file
  buffer <- matrix (scan (files [1], ...), ncol = 2, byrow = TRUE)</pre>
  ## first column gives the wavelength vector
 wavelength <- buffer [, 1]</pre>
  ## preallocate the spectra matrix:
 ## one row per file x as many columns as the first file has
  spc <- matrix (ncol = nrow (buffer), nrow = length (files))</pre>
  ## the first file's data goes into the first row
  spc [1, ] <- buffer [, 2]</pre>
  ## now read the remaining files
  for (f in seq (along = files)[-1]) {
   buffer <- matrix (scan (files [f], ...), ncol = 2, byrow = TRUE)
```

```
## check whether they have the same wavelength axis
if (! all.equal (buffer [, 1], wavelength))
    stop (paste(files [f], "has different wavelength axis."))

spc [f, ] <- buffer[, 2]
}

## make the hyperSpec object
new ("hyperSpec", wavelength = wavelength, spc = spc,
    data = data.frame (file = files), label = label)
}</pre>
```

Note how the labels are set. The label with the special name .wavelength corresponds to the wavelength axis, all data columns should have a label with the same name. The spectra are always in a data column called spc.

Thus,

imports the spectra.

```
> source ("scan.txt.PerkinElmer.R")
> scan.txt.PerkinElmer ("txt.PerkinElmer/flu?.txt", skip = 54)
hyperSpec object
   6 spectra
   2 data columns
   181 data points / spectrum
wavelength: lambda/nm [numeric] 405.0 405.5 ... 495
data: (6 rows x 2 columns)
   1. file: [factor] txt.PerkinElmer/flu1.txt txt.PerkinElmer/flu2.txt ... txt.PerkinElmer/flu6.txt
   2. spc: I[f1]/"a.u." [matrix181] 27.150 66.801 ... 294.65
```

This function is not exported by *hyperSpec*. While it is already useful for importing files, it is not yet general enough to work immediately with new data: the file header is completely ignored. Thus information like the excitation wavelength is lost.

# 7.2. Deriving a More Specific Function: read.ENVI.Nicolet

The function read. ENVI.Nicolet is a good example for a more specific import filter derived from a general filter for the respective file type. Nicolet FT-IR Imaging software saves some non-standard keywords in the header file of the ENVI data. These information can be used to reconstruct the x and y axes of the images. The units of the spectra are saved as well.

read.ENVI.Nicolet thus first adjusts the parameters for read.ENVI. Then read.ENVI does the main work of importing the file. The resulting *hyperSpec* object is post-processed according to the special header entries.

For using the function, see section 6.3 (p. 14).

```
read.ENVI.Nicolet.R

read.ENVI.Nicolet <- function (..., # goes to read.ENVI

# file headerfile, header

x = NA, y = NA, # NA means: use the specifications from the header file if possible

log = list (),

keys.hdr2log = TRUE,

nicolet.correction = FALSE) {

## set some defaults

log <- modifyList (list (short = "read.ENVI.Nicolet",

long = list (call = match.call ())),
```

```
log)
   ## the additional keywords to interprete must be read
   if (! isTRUE (keys.hdr2log))
       keys.hdr2log <- unique (c ("description", "z plot titles", "pixel size", keys.hdr2log))
   ## most work is done by read.ENVI
   spc <- read.ENVI (..., keys.hdr2log = keys.hdr2log,</pre>
                                       x = if (is.na(x)) 0 : 1 else x,
                                       y = if (is.na (y)) 0 : 1 else y,
                                       log = log)
   ## get the header for post-processing
   header <-spc@log$long.description [[1]]$header
### From here on processing the additional keywords in Nicolet's ENVI header *********************
   ## z plot titles ------
   ## default labels
   label <- list (x = expression (\dot{\ }/\dot{\ } (x, micro * m)),
                                 y = expression ('/' (y, micro * m)),
                                 spc = 'I / a.u.',
                                  .wavelength = expression (tilde (nu) / cm^-1))
   ## get labels from header information
   if (!is.null (header$'z plot titles')){
       pattern <- "^[[:blank:]]*([[:print:]^,]+)[[:blank:]]*,.*$"</pre>
        tmp <- sub (pattern, "\\1", header$'z plot titles')</pre>
       if (grepl ("Wavenumbers (cm-1)", tmp, ignore.case = TRUE))
           label$.wavelength <- expression (tilde (nu) / cm^(-1))</pre>
           label $. wavelength <- tmp
       pattern <- "^[[:blank:]]*[[:print:]^,]+,[[:blank:]]*([[:print:]^,]+).*$"</pre>
       tmp <- sub (pattern, "\\1", header$'z plot titles')
if (grepl ("Unknown", tmp, ignore.case = TRUE))</pre>
          label$spc <- "I / a.u."
        else
           label$spc <- tmp</pre>
   ## modify the labels accordingly
   spc@label <- modifyList (label, spc@label)</pre>
   ## set up spatial coordinates -----
   ## look for x and y in the header only if x and y are NULL
   ## they are in `description` and `pixel size`
   ## set up regular expressions to extract the values
   {\tt p.description} \begin{tabular}{ll} {\tt p.description} &\leftarrow {\tt paste} & ("`Spectrum position [[:digit:]]+ of [[:digit:]]+ positions,", the property of the pr
                                                   "X = ([[:digit:].-]+), Y = ([[:digit:].-]+)$")
   p.pixel.size <- "^[[:blank:]]*([[:digit:].-]+),[[:blank:]]*([[:digit:].-]+).*$"
    if (is.na (x) && is.na (y) &&
           ! is.null (header$description) && grepl (p.description, header$description ) &&
           ! is.null (header$'pixel size') && grepl (p.pixel.size, header$'pixel size')) {
       x [1] <- as.numeric (sub (p.description, "\\1", header$description)) y [1] <- as.numeric (sub (p.description, "\\2", header$description))
       x [2] <- as.numeric (sub (p.pixel.size, "\\1", header^ppixel size')) y [2] <- as.numeric (sub (p.pixel.size, "\\2", header^ppixel size'))
       \hbox{\it \#\# it seems that the step size is given in } \hbox{\it mm while the offset is in micron}
       if (nicolet.correction) {
```

```
x [2] <- x [2] * 1000
y [2] <- y [2] * 1000
}

## now calculate and set the x and y coordinates
x <- x [2] * spc$x + x [1]
if (! any (is.na (x)))
    spc@data$x <- x

y <- y [2] * spc$y + y [1]
if (! any (is.na (y)))
    spc@data$y <- y
}
spc@data$y <- y
}</pre>
```

# 7.3. Deriving import filters for spc files

Please note that future changes inside the read.spc function are likely to occur. However, if you just post-process the *hyperSpec* object returned by read.spc, you should be fine.

#### References

- [1] Robert S. McDonald and Jr. Paul A. Wilks. Jcamp-dx: A standard form for the exchange of infrared spectra in computer readable form. *Applied Spectroscopy*, 42(1):151–162, 1988.
- [2] Henrik Bengtsson. R.matlab: Read and write of MAT files together with R-to-MATLAB connectivity, 2014. URL http://CRAN.R-project.org/package=R.matlab. R package version 3.1.1.
- [3] Universal Data Format Specification. Galactic Industries Corp., 1997. URL http://ftirsearch.com/features/converters/gspc\_udf.zip.

# A. File Import Functions by Format

Format	Manufacturer	Function	section	Notes
$ASCII\ file\ formats$				
ASCII long		read.txt.long	4, p. 5	
ASCII long	Renishaw (Raman)	scan.txt.Renishaw	6.6.1, p. 16	
ASCII long	Kaiser (Raman)	read.txt.long	6.5.1, p. 15	Not recommended, see discussion
ASCII long	Perkin Elmer (Fluorescence)	read.txt.PerkinElmer	7.1, p. 19	Reads multiple files, needs to be sourced.
ASCII wide		read.txt.wide	4, p. 5	
ASCII wide	Horiba Jobin Yvon	read.txt.wide	6.7, p. 17	e. g. LabRAM spectrometers
ASCII wide transposed	Witec (Raman)	scan.txt.Witec	6.8, p. 18	Export Table
JCAMP-DX		-	4.2, p. 6	see read.jdx
JCAMP-DX	Renishaw (Raman)	-	4.2, p. 6	not available
JCAMP-DX	Shimadzu (GC,GC-MS)	jcamp-dx	4.2, p. 6	import for subset of the JCAMP-DX standard
JCAMP-DX	PerkinElmer (Infrared)	jcamp-dx	4.2, p. 6	import for subset of the JCAMP-DX standard
Witec ASCII	Witec (Raman)	scan.dat.Witec	6.8, p. 18	Save ASCII X, Save ASCII Y
binary file formats				
array		-	2, p. 3	
ENVI		read.ENVI	5.2, p. 10	
ENVI	Bruker (Infrared Imaging)	read.ENVI	6.2, p. 14	
ENVI	Nicolet (Infrared Imaging)	read.ENVI.Nicolet	6.3, p. 14	
hol	Kaiser (Raman)	-	6.5, p. 15	via Matlab
Matlab	Matlab	R.matlab::readMat	5.1, p. 8	
Matlab	Cytospec	read.cytomat	5.1.2, p. 9	
matrix		-	2, p. 3	
Opus	Bruker (Infrared Imaging)	-	6.2, p. 14	
spc		read.spc	5.3, p. 10	
spc	Kaiser (Raman Map)	read.spc.KaiserMap	6.5.2, p. 15	Reads multiple files
$\operatorname{spc}$	Kaiser (Raman)	read.spc	5.3, p. 10	Reads multiple files
$\operatorname{spc}$	Renishaw (Raman)	read.spc	6.6.1, p. 16	Not recommended, see discussion of ASCII files.
spc	Witec (Raman)	read.spc	5.3, p. 10	spc export not available for images

# B. File Import Functions by Manufacturer

Manufacturer	Format	Function	section	Notes
Manufacturers				
Bruker (Infrared Imaging)	ENVI	read.ENVI	6.2, p. 14	
Bruker (Infrared Imaging)	Opus	-	6.2, p. 14	
Cytospec	Matlab	read.cytomat	5.1.2, p. 9	
Horiba Jobin Yvon	ASCII wide	read.txt.wide	6.7, p. 17	e. g. LabRAM spectrometers
Kaiser (Raman)	ASCII long	read.txt.long	6.5.1, p. 15	Not recommended, see discussion
Kaiser (Raman)	hol	-	6.5, p. 15	via Matlab
Kaiser (Raman Map)	$\operatorname{spc}$	read.spc.KaiserMap	6.5.2, p. 15	Reads multiple files
Kaiser (Raman)	$\operatorname{spc}$	read.spc	5.3, p. 10	Reads multiple files
Matlab	Matlab	R.matlab::readMat	5.1, p. 8	
Nicolet (Infrared Imaging)	ENVI	read.ENVI.Nicolet	6.3, p. 14	
PerkinElmer (Infrared)	JCAMP-DX	jcamp-dx	4.2, p. 6	import for subset of the JCAMP-DX standard
Perkin Elmer (Fluorescence)	ASCII long	read.txt.PerkinElmer	7.1, p. 19	Reads multiple files, needs to be sourced.
Renishaw (Raman)	ASCII long	scan.txt.Renishaw	6.6.1, p. 16	
Renishaw (Raman)	JCAMP-DX	-	4.2, p. 6	not available
Renishaw (Raman)	$\operatorname{spc}$	read.spc	6.6.1, p. 16	Not recommended, see discussion of ASCII files.
Shimadzu (GC,GC-MS)	JCAMP-DX	jcamp-dx	4.2, p. 6	import for subset of the JCAMP-DX standard
Witec (Raman)	Witec ASCII	scan.dat.Witec	6.8, p. 18	Save ASCII X, Save ASCII Y
Witec (Raman)	ASCII wide transposed	scan.txt.Witec	6.8, p. 18	Export Table
Witec (Raman)	spc	read.spc	5.3, p. 10	spc export not available for images

# C. File Import Functions by Spectroscopy

Spectroscopy	Format	Manufacturer	Function	section	Notes
Fluorescence	ASCII long	Perkin Elmer	read.txt.PerkinElmer	7.1, p. 19	Reads multiple files, needs to be sourced.
GC,GC-MS	JCAMP-DX	Shimadzu	jcamp-dx	4.2, p. 6	import for subset of the JCAMP-DX standard
Infrared	JCAMP-DX	PerkinElmer	jcamp-dx	4.2, p. 6	import for subset of the JCAMP-DX standard
Infrared Imaging	ENVI	Bruker	read.ENVI	6.2, p. 14	
Infrared Imaging	ENVI	Nicolet	read.ENVI.Nicolet	6.3, p. 14	
Infrared Imaging	Opus	Bruker	-	6.2, p. 14	
Raman	ASCII long	Renishaw	scan.txt.Renishaw	6.6.1, p. 16	
Raman	ASCII long	Kaiser	read.txt.long	6.5.1, p. 15	Not recommended, see discussion
Raman	ASCII wide transposed	Witec	scan.txt.Witec	6.8, p. 18	Export Table
Raman	hol	Kaiser	-	6.5, p. 15	via Matlab
Raman	JCAMP-DX	Renishaw	-	4.2, p. 6	not available
Raman	$\operatorname{spc}$	Kaiser	read.spc	5.3, p. 10	Reads multiple files
Raman	spc	Renishaw	read.spc	6.6.1, p. 16	Not recommended, see discussion of ASCII files.
Raman	spc	Witec	read.spc	5.3, p. 10	spc export not available for images
Raman	Witec ASCII	Witec	scan.dat.Witec	6.8, p. 18	Save ASCII X, Save ASCII Y
Raman Map	$\operatorname{spc}$	Kaiser	read.spc.KaiserMap	6.5.2, p. 15	Reads multiple files

# Index

Agilent	Map, 20
ENVI, 15	Nicolet, 20
ASCII	PerkinElmer, 6
compressed, 16	initialize hyperSpec object, 3
JCAMP-DX, 6	
long, 5	JCAMP-DX
Fluorescence, 19	ASCII, 6
Kaiser, 15	Infrared, 6
PerkinElmer, 19	PerkinElmer, 6
Raman, 15, 16	Shimadzu, 6
Renishaw, 16	jdx, see JCAMP-DX
samples in columns, 5	
transposed, 5	Kaiser
wide, 5	ASCII long, 15
Horiba, 17	hol, 15
Horiba Jobin Yvon, 17	Map, 15
LabRAM, see Horiba	spc, 10, 15
Witec, 18	• , ,
zip, 16	LabRAM, see Horiba
atomic emission	Edistriivi, occ Horista
NIST, 8	Mon. 15
11101,0	Map, 15
Bruker	ENVI, 10, 14, 15
	Kaiser, 15
AXS, 5	Raman, 15
ENVI, 10, 14	Matlab, 8
powder diffraction, 5	Cytospec, 9
x-ray, 5	
	new hyperSpec object, 3
create hyperSpec object, 3	Nicolet
Cytospec, 9	ENVI, 14, 20
Matlab, 9	Infrared, 20
	Map, 20
ENVI	NIST
Agilent, 15	atomic emission, 8
Bruker, 10, 14	atomic chibbion, o
Infrared, 10, 14, 15	Dl.: El
Nicolet, 20	PerkinElmer
Map, 10, 14, 15, 20	ASCII long, 19
	Fluorescence, 19
Nicolet, 14, 20	Infrared, 6
Varian, 10	JCAMP-DX, 6
	powder diffraction
Fluorescence	Bruker, 5
ASCII long, 19	
PerkinElmer	Raman
ASCII, 19	ASCII long, 15, 16
FT-IR, see Infrared	hol, 15
	Horiba
hol	ASCII wide, 17
Kaiser, 15	Horiba Jobin Yvon
Horiba	ASCII wide, 17
ASCII	Kaiser, 10, 15
wide, 17	LabRAM, see Horiba
Horiba Jobin Yvon	Map, 15
ASCII	Renishaw, 10
wide, 17	ASCII, 16
hyperSpec object	
create, 3	spc, 16
Crowner, o	spc, 10, 15
Imaga aga Man	Witec
Image, see Map	ASCII, 18
Infrared	Export Table, 18
ENVI, 10, 14, 15, 20	Save ASCII X, Save ASCII Y, 18
JCAMP-DX, 6	spc, 18

```
reference
NIST
atomic emission, 8
Renishaw
ASCII long, 16
Raman, 16
spc, 10, 16

Shimadzu
GC-MS, 6
JCAMP-DX, 6
spc, 10
Kaiser, 10, 15
Raman, 10, 15, 16
Renishaw, 10, 16
TriVista, 10
Witec, 18

TriVista
spc, 10

Varian, see Agilent
ENVI, 10

Witec
ASCII, 18

x-ray
Bruker, 5
```

# **Session Info**

R version 3.1.2 (2014-10-31)

Platform: x86\_64-pc-linux-gnu (64-bit)

locale:

[7] LC\_PAPER=de\_DE.UTF-8 LC\_NAME=C LC\_ADDRESS=C

[10] LC\_TELEPHONE=C LC\_MEASUREMENT=de\_DE.UTF-8 LC\_IDENTIFICATION=C

attached base packages:

[1] grid stats graphics grDevices utils datasets methods base

other attached packages:

[1] R.matlab\_3.1.1 hyperSpec\_0.98-20150217 mvtnorm\_1.0-2 ggplot2\_1.0.0

[5] lattice\_0.20-30

loaded via a namespace (and not attached):

[1] colorspace\_1.2-4 digest\_0.6.8 gtable\_0.1.2 MASS\_7.3-39 munsell\_0.4.2 [6] plyr\_1.8.1 proto\_0.3-10 Rcpp\_0.11.4 reshape2\_1.4.1 R.methodsS3\_1.6.1 [11] R.oo\_1.18.0 R.utils\_1.34.0 scales\_0.2.4 stringr\_0.6.2 tools\_3.1.2